

INTENDED USE


Product for use in the quantitative determination in vitro of the γ -GT activity in human serum. The results of the test must always be interpreted in conjunction with the clinical picture.

INTRODUCTION

Although γ -GT IFCC is present in a large number of tissues, the enzyme which we wish to detect in the serum is principally part of the hepato-biliary system. Consequently raised levels of γ -GT are seen in all forms of disease or damage to the hepatic system. From the clinical viewpoint, the enzyme is useful for the diagnosis of obstructive jaundice, cholangitis and cholecystitis. High γ -GT levels are also seen during the assumption of alcohol or pharmaceutical products (sedatives, anticonvulsants and tranquilizers) (1).

PRINCIPLE

Gamma glutamyl transpeptidase (y-GT) catalyzes the transfer of the g-glutamyl group from the substrate g-glutamyl-3-carboxy-4-nitroanilide to glycylglycine releasing L-glutamyl-glycylglycine and 5-amino-2-nitrobenzoate.




The rate of formation of 5-amino-2-nitrobenzoate, determined kinetically at 405 nm, is proportional to y-GT activity.

REAGENTS

- Kit R1: 3 x 20 mL R2: 3 x 20 mL code A-R0000000013
Reagent 1: no. 3 vials x 19.5 mL code 184025
Reagent 2: no. 3 vials x 5.0 mL code 184026
 Kit R1: 4 x 50 mL R2: 4 x 20 mL code A-R0200000601
Reagent 1: no. 4 vials x 40.5 mL code 184052
Reagent 2: no. 4 vials x 10.5 mL code 184053

Concentrations

Reagent 1:			
	Conc.	U.M.	
y-glutamyl-3-carboxi-4-nitroanilide	4.00	mM	
Reagent 2:			
Glycylglycine	750	mM	

Precautions

Reagent B is harmful and must be used with caution. Refer to the appropriate Material Safety Data Sheet for the risk phrases and safety measures (**Identification of Risks: GHS07 - H315, H319**). In addition to any risk indications related to the active components, reagents may contain inactive components such as preservatives and detergents. The total concentrations of these components are lower than the limit reported in the EC directive 1272/2008 and subsequent amendments and additions. However, It is recommended to handle reagents according to the standards of good laboratory practice and adopt all adequate protective measures.

Storage and stability

Store at 2 - 8°C and protect from direct light. When correctly stored, the reagents are stable up the expiry date reported on the label. After opening, the vial R1 and R2 are stable 30 days if recapped immediately and protected from contamination, evaporation, direct light and stored at correct temperature.

SAMPLE COLLECTION

Type of sample and storage

Use serum or plasma with EDTA, heparin, citrate or oxalate/fluoride (5). Use serum free of hemolysis. y-GT is stable in serum for at least a week if stored at - 4°C to + 20°C (6) and 3 months at - 20°C (5).

Precaution

All human samples must be handled and disposed of as potentially infectious materials.

Procedure

Quality control

Human control serum with known levels of y-GT is commercially available for quality control purposes. Data sheets are included, listing the values and the confidence limits. Normal and abnormal control sera are available from I.S.E. S.r.l. "Normal control serum" code R0400000006 and "Pathological control serum" code R0400000106. Obtained values must be within the range of acceptability. If erratic results occur, the following points should be checked:

- Cleanliness of glassware.
- Wavelength setting.
- Expiration date of reagents.

Automation

This kit, though developed and manufactured to be used as manual assay and with I.S.E. S.r.l. analyzer, can be used also with other analyzers able to meet the specifications indicated in section "Reaction conditions - Test procedure" Application sheets are available for automatic instruments.

All applications not explicitly approved by I.S.E. S.r.l. cannot be guaranteed in terms of performance, and must therefore be established by the operator.

Calibration

For calibration use the "Multicalibrator" I.S.E. S.r.l. code R0300000006.

Traceability:

The y-GT concentration is reported in the package insert supplied with the Calibrator Serum.

Calibration Stability

For the instrumentation series Miura, the calibration is recommended to be done every 10 days.

Method for automated instrumentation

Analyzer:	I.S.E. Miura		
Analyte Name :	Gamma GT	Ref.:	A-R0000000013 A-R0200000601
Analyte Code:	γ GT		
Method:	Kinetic		
Unit:	U/L		
Filter F1:	405 nm		
Blank in:	Not Used		
Step	Reaction volume	U.M.	
Volume reagent R1:	160	μ L	
Volume reagent R2:	40	μ L	
Sample volume:	20	μ L	
Final Incubation:	60	Sec.	
Kinetic reading time	240	Sec.	

Materials included in the kit

The reagents are described above.

Materials required but not supplied in the kit

Calibrators and controls.

NORMAL VALUE

Serum:

- Male: 11 – 51 U/L
- Female: 7 – 33 U/L

Each laboratory should calculate its own normal values on the basis of its local population

ANALYTICAL CHARACTERISTICS / PERFORMANCE

Linearity

If concentration at 405 nm is greater than > 500 IU/L, repeat the test using serum diluted 1:10 with isotonic saline. Multiply by 10.



Specificity

Under the conditions of the assay system this method is specific for γ -glutamyl transpeptidase.

Accuracy-Recovery

The recovery of pure γ -GT added to normal sample at known titer was 101.7%.

Precision of the method

Within-run precision					
Range	U.M.	Mean	S.D.	C.V. (%)	No.
Low	U/L	36.53	1.311	3.59%	18
High	U/L	152.7	3.774	2.47%	18
Between-run precision					
Range	U.M.	Mean	S.D.	C.V. (%)	No.
Low	U/L	36.53	0.756	2.07%	18
High	U/L	152.7	1.317	0.86%	18

Sensitivity

At λ 405 nm a concentration of about 9,00 U/L of γ -GT in the conditions established for this test.

Comparative method

The I.S.E. S.r.l. method was compared to a method in use. A linear regression slope $Y = 1.017x + 5.09$ was determined. A Correlation coefficient $r = 1.000$ was calculated. Sample tested No. = 48.

Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.


The product is in conformity with D.L. 8 September 2000, no. 332 "Actuation of the directive 98/79/EC regarding in vitro medical diagnostic devices".

Symbols on labels and packaging


 = In vitro diagnostic medical device


 = Catalog Number

 = Lot Number

 = Manufacturer

 = Expiration date

 = Temperature limitation

 = Instruction for use

References

- Orlowski M., Meister A. Biochim Biophys Acta; 73: 679,1963.
- Guder W G, Narayanan s, wiser H, Zawta b. List of Analytes; Preanalytical variables. Brochure in: Samples: From Patient to the Laboratory. GIT Verlag GmbH, Darmstadt, 1996.
- IFCC 2002/7: IFCC Primary Reference Procedures for the Measurement of Catalytic Activity Concentrations of Enzymes at 37°C. Part. 6. Reference Procedure for the Measurement of Catalytic Concentration of γ -Glutamyltransferase, Clin. Chem. Lab. Med. 2002, 40(7): 734-738.
- Nielsen L.G. et al. Amer J Med Technol; 44: 279,1978

