LOT Lot:	Exp.:	Storage: 2-8 °C
-----------------	-------	-----------------

Intended use

For calibration of the I.S.E. S.r.l. test methods on all the I.S.E. S.r.l. Clinical Chemistry Analyzers. The reagents and methods described are manufactured and controlled by "I.S.E. S.r.l."

REAGENT INFORMATION

Kit 3 x 1 mL Code R030000009
Lyophilized 1: N° 3 x 1 mL vials
Nature of the product lyophilized, obtained from a pool of human sera code 184253

Each calibrator serum (1x1mL) contains HDL and LDL cholesterol by human matrix. The concentration of CHO-HDL/LDL is obtained by the direct method for selective inhibition. (Reag. ISE Srl, CHO-HDL code. A-R0200001301; CHO- LDL cod. A-R0200001401)

Analyte	Assigned Value	Units
HDL	58,3	mg/dL
	1,51	mmol/L
LDL	154,0	mg/dL
	3,99	mmol/L

Stability

Calibrator serum, stored at 2-8 °C in a sealed vial is stable until the expiration date printed on the package.

The reconstituted material is stable for 4 days at 2-8 °C or 1 month at -20 °C (If thawed only once).

Method of use

- Remove the serum from the refrigerator.
- Remove the cap and add exactly 1.0 ml of distilled water. Replace the cap, then mix gently by inversion to avoid foaming.
- Let stand 10 minutes before performing the calibration of HDL Cholesterol
- Let stand 2 hours before performing the calibration of LDL Cholesterol
- Close the serum and store at 2-8 °C or -20 °C taking into account the stability attributed. After an initial thaw do not refreeze the product.

Warnings

Occasionally, the serum may have slightly cloudy. This is a feature of the lyophilized material. We recommend the use as specified in the recommended applications.

Precautions and warnings

The sera have been tested and proved negative for the presence of HbsAg, and also for anti-HIV 1 and 2 and anti-HCV antibodies. However, the material must be considered potentially dangerous and handled in the same way as the test samples.

Disposal

Disposal must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

The product is in conformity with D.L. 8 September 2000, no. 332 "Actuation of the directive 98/79/CE regarding in vitro medical diagnostic devices".

Symbols used on labels and packaging

IVD = In vitro diagnostic medical device

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I.S.E S.r.l.