

6 x 50 mL • REF A-R0100000301

3 x 18,5 mL • REF R3330000012

INTENDED USE

Product for use in the quantitative determination in vitro of the Calcium concentration in human serum, plasma or urine. The results of the test must always be interpreted in conjunction with the clinical context. For professional use only.

INTRODUCTION

Calcium is one of the mineral elements present in the blood, and one of the chief constituents of bone tissue. It is present in small amounts in soft tissues and in extracellular fluids. About 50% of the calcium present in the blood is free, 40% is bound to proteins and 10% is present in compounds. About 80% of calcium bound to proteins is associated with albumin and the remaining 20% with globulins. Extracellular calcium has the role of maintaining intracellular calcium, mineralizing bone, blood coagulation and maintenance of the cell membrane. Calcium is of particular importance in muscular contraction, in addition to being a second messenger which influences enzymatic activity and hormonal secretion. High serum calcium levels (hypercalcemia) are essentially due to an increased flow of extracellular calcium produced by bone, intestinal and renal tissue. Low serum calcium levels (hypocalcemia) may be due to a reduction in the calcium bound to albumin or to the free fraction.

PRINCIPLE

Calcium ions form a highly coloured complex with Arsenazo III at neutral pH. The colour which has formed will correspond to the amount of calcium present in the sample.

Calcium + 2 Arsenazo III → Ca-Arsenazo Complex

REAGENTS


A-R0100000301 – 6 x 50 mL

Reagent: no. 6 vials x 50.0 mL ready to use

R3330000012 – 3 x 18.5 mL

Reagent: no. 3 vials x 14.0 mL ready to use

Concentrations

| Reagent: | Conc. | U.M. |  GHS08* |
|---------------------------|-------|------|--|
| Imidazol buffer (pH 6.75) | 100 | mM | |
| Arsenazo III | 0.20 | mM | |
| Stabilizing agents | - | - | |

* Warning **DANGER** Contains: imidazole (CAS 288-32-4)

H360D- May damage the unborn child.

P201 - Obtain special instructions before use.

P202 - Do not handle until all safety precautions have been read and understood.

P280 - Wear protective gloves/protective clothing/eye protection/face protection/hearing protection.

P308+P313 - IF exposed or concerned: Get medical advice/attention.

P501 - Dispose of contents/container in accordance with local/ Regional /national /international regulations.

Precautions

Kit for professional laboratory use, used only by qualified and properly trained technical personnel, under the supervision of a doctor in charge of the laboratory. In addition to risk indications related to the active components, reagents may contain inactive components such as preservatives and detergents. The total concentration of these components is lower than the limits reported in the EC 1272/2008 Regulation and subsequent amendments and additions. However, it is recommended to handle reagents according to the standards of good laboratory practice

Reports of serious incidents

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device.

Storage and stability

Store at 2 - 8°C and protect from direct light. When correctly stored, the reagents are stable up to the expiry date reported on the label. A slight variation in the composition of the reagents may occur between batches, but this has no effect on the test results.

After opening, the vial R1 and R2 are stable 30 days if recapped immediately and protected from contamination, evaporation, direct light and stored at correct temperature.

Reagent Preparation

The reagent is liquid ready for use. After opening, the reagent is stable for 30 days if closed and stored at 2 - 8°C. Do not mix different batches.

SAMPLE COLLECTION

Type of sample and storage

Heparinized plasma or serum samples should be used. Do not use anticoagulants containing citrates, oxalates or EDTA as these tend to remove the calcium through formation of complexes. The samples can be stored for 7 days at 20 - 25°C, 3 weeks at 4 - 8°C, 8 months at - 20°C.

Precautions

All human samples must be handled and disposed of as potentially infectious materials.

PROCEDURE

Quality control

Control sera with known Calcium content are commercially available for quality control, including certificates of analysis showing the values and limits of confidence. Normal and pathological control sera are available as "Normal control serum" cod. R0400000006 and "Pathological control serum" code R0400000106. The values obtained must be contained within the acceptability range. In case of incorrect results check the following points:

- Cleanliness of glassware.
- Wavelength setting.
- Expiration date of reagents.

Automation

Although this device has been developed and manufactured to be used by manual method and instrumental systems, it can also be used in combination with other instrumental devices that are able to meet the requirements indicated in the paragraph "Reaction conditions / Technique".

All applications that are not explicitly approved cannot be guaranteed in terms of performance and will therefore have to be evaluated by the user.

Calibration

For calibration use the kit "Multicalibrator" code R0300000006. For ISE srl instruments, calibration is recommended every 10 days.

Traceability

Calcium traceability is visible in the insert of the calibration serum packaging.

MANUAL METHOD

The kit, in Open format, can be used manually through the use of spectrophotometer or photometer with the following parameters:

Reaction conditions

Wavelength (primary): 660 nm
Temperature: 37°C

Technical – procedure with Serum as starter

Bring the reagents to reaction temperature and operate away from direct light.

| | U.M. | Blank | Calibr. Serum | Sample |
|---------------|------|-------|---------------|--------|
| Reagent | µL | 1000 | 1000 | 1000 |
| Calibr. serum | µL | - | 20 | - |
| Sample | µL | - | - | 20 |
| Blank | µL | 20 | - | - |

Mix well and incubate for 5 minutes at 37°C.

Measure absorbance of the sample and standard against reagent blank.

The reaction volumes may be varied proportionally without alteration of results.

Results:

Manual Method

Calculation of Calcium:

$$\frac{\text{O.D. Sample}}{\text{O.D. Calibrator serum}} \times \text{Calibrator serum Concentration} = \text{Calcium mg/dL}$$

NORMAL VALUES

Serum or plasma: 8.6 – 10.3 mg/dL (2.15 – 2.57 mmol/L)

Urine: 100 – 400 mg/24 h (2.5 – 10 mmol/24 h)

Each laboratory must establish its own normal values on the basis of its local catchment area.



Calcium All – Instructions for use (IFU)

Distributed by:

ISE S.r.l.

CUSTOMISED SOLUTIONS
FOR YOUR LABORATORY

IVD

CE

6 x 50 mL •  A-R0100000301

3 x 18,5 mL •  R3330000012

Reagents included in the kit

The reagents are described above.

Materials required but not supplied in the kit

Calibrators and controls.

ANALYTICAL CHARACTERISTICS/PERFORMANCE

Linearity

The method is linear up to 17.90 mg/dL.

Specificity

The method is specific for the determination of Calcium in the analytical conditions described.

Accuracy – Recovery

The recovery of Calcium added to normal samples at known concentrations showed an accuracy of 94.4%.

Interferences

The high dilution of the sample with the reagent reduces to a minimum the interference by lipids. Bilirubin below 58 mg/dL does not interfere in the reaction, hemoglobin interferes at concentrations above 10 g/L. Copper may cause interference. It is therefore advisable to avoid contamination with reagents containing high concentrations of copper such as Biuret for the reading of total proteins. For other interfering substances, make reference to the bibliography reported below (4).

Precision of the method

| Within-run precision | | | | | |
|-----------------------|-------|-------|-------|----------|-----|
| Range | U.M. | Mean | S.D. | C.V. (%) | No. |
| Low | mg/dL | 5.46 | 0.166 | 3.03% | 20 |
| High | mg/dL | 11.65 | 0.133 | 1.14% | 20 |
| Between-run precision | | | | | |
| Range | U.M. | Mean | S.D. | C.V. (%) | No. |
| Low | mg/dL | 5.46 | 0.081 | 1.49% | 20 |
| High | mg/dL | 11.65 | 0.754 | 6.47% | 20 |

Sensitivity

At λ 630 nm, a concentration of about 0.1 mg/dL of Calcium in the conditions established for this test.

Comparative method

The Calcium Arsenazo III reagent was compared with a similar commercial method. Samples tested No. = 71; $Y = 0.9069x + 0.97$; $r = 0.989$.

Disposal of reagents

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Manufacturer:

Scavo Diagnostics International

Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)

Phone +39 0577 39041 - Fax +39 0577 390 444

Distributor:

I.S.E S.r.l.








Via Delle Driadi, 45 – 00133 Roma

Tel.+39 077 4579365; FAX +39 077 4579305

E-mail: info@logotech-ise.com

www.logotech-ise.com

Symbols used in IFU and Packaging

| | |
|--|--|
|  In vitro diagnostic medical device vitro |  Manufacturer |
|  Catalogue Number |  Instruction for use |
|  Lot Number |  Temperature limitation |
|  Expiration date | |

Reference

1. Butris CA and Ashwood ER (Ed.), Tietz Fundamentals of Clinical Chemistry, 5th Edition, W.B. Saunders Company, Philadelphia, 2001, p.797-799, 968.
2. Janssen JW and Helbing AR. Arsenazo III. An improvement of the routine calcium determination in serum. Eur. J. Clin. Chem. Clin. Biochem. 29 (3) pp. 197-201, 1991.
3. Guder WG, Narayanan S., Wisser H., Zavata B. List of analyses; preanalytical variables. Brochure in: Samples: from patient to the laboratory. Git Verlag GmbH, Darmstadt, 1996.
4. Young D, Effects of drugs on clinical laboratory tests. 5th Edition, AACC Press, Washington, DC, 3-149 – 3-158, 2000.

| REVISION | DATE | CHANGE |
|----------|---------|--|
| Rev.A | 12/2022 | New Issue for IVDR Regulation (UE) 2017/746 compliance |



Scavo Diagnostics International
Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)
Phone +39 0577 39041 - Fax +39 0577 390 444