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**ISE S.r.l.**  
CUSTOMISED SOLUTIONS  
FOR YOUR LABORATORY

R1: 4 x 5 mL •  R040000006

### Intended use

The Normal Serum is used as a quality control reagent to test the accuracy and precision of quantitative diagnostic methods on I.S.E. S.r.l. biochemical analyzers. The theoretical control values are listed on the product information sheet enclosed. The concentrations or theoretical activities described are given as an average value within the normal range. For professional use only.

### Composition

The Normal Serum is obtained from defibrinated, freeze-dried human plasma. The declared nature of the product denotes the potential risk involved and it must therefore be used, adopting the same precautions used in handling the test samples.

### Reconstitution of the lyophilized serum

Remove the cap with care, to avoid losing traces of the freeze-dried product. Add exactly 5 mL of distilled water, recap, mix by inversion a few times and leave to rest for about 30 minutes. Before use, mix again by inversion, avoiding foam formation. **Important:** when testing for Alkaline Phosphatase (ALP), wait a further 60 minutes until the enzyme is completely reactivated.

### Storage and stability of the lyophilized serum

Store at 2 - 8°C. The stability of the serum before reconstitution is reported on the vial label.

Stability of the reconstituted serum:

at 18 - 25°C 8 hours  
at 2 - 8°C 1 week, 2 days for Bilirubin and 1 day for Acid Phosphatase  
at (-15) - (-25)°C 1 month (If thawed only once)

### Origin of the biological components

Analyte	Origin
ALT (GPT)	Pig heart
AST (GOT)	Pig heart
Albumin	Human serum
Amylase	Human saliva/pig pancreas
Cholesterol	Bovine plasma
Cholinesterase	Human serum
CK	Pig plasma
Y-GT	Pig kidney
LDH	Pig heart
ALP (Alkaline phosphatase)	Calf intestine
ACP (acid phosphatase)	Human prostate/potato
Total Proteins	Human serum

### Precautions and warnings

The sera have been tested and proved negative for the presence of HbsAg, and also for anti-HIV 1 and 2 and anti-HCV antibodies. However, the material must be considered potentially dangerous and handled in the same way as the test samples.

### Reporting serious incidents

Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union where the user and/or patient is established, of cases of serious accidents that have occurred in relation to the device. For other jurisdictions, reports of serious incidents must be made in accordance with regulatory requirements. By reporting serious incidents, you help provide more information about the safety of the in vitro diagnostic medical device.

### Assigned values and traceability

The theoretical values have been calculated on the basis of the results obtained on I.S.E. S.r.l. biochemical analyzers using I.S.E. S.r.l. reagents. The values for each single method have been obtained using the I.S.E. S.r.l. Multicalibrator code R030000006. The acceptable range with reference to the theoretical value was calculated within  $\pm 3$  standard deviations. The traceability of the values, with reference to each diagnostic method, is reported on the enclosed product information sheet. The use of procedures or reagents other than those indicated may lead to differences between the expected values and those obtained.

### Materials included in the kit

- Freeze-dried serum as described
- Product information sheet reporting the theoretical values of the control serum

### Materials required but not included in the kit

- Distilled or deionized water
- Reagents for I.S.E. S.r.l. analytical systems

### Disposal

Disposal must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

### Manufacturer:

**Scavo Diagnostics International**  
Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)  
Phone +39 0577 39041 - Fax +39 0577 390 444

### Distributor:

**I.S.E S.r.l.**  
Via Delle Driadi, 45 – 00133 Roma  
Tel. +39 077 4579365; FAX +39 077 4579305  
E-mail: info@logotech-ise.com  
www.logotech-ise.com

### Symbols used on labels and packaging

 = In vitro diagnostic medical device

 = Catalogue Number

 = Lot Number

 = Manufacturer

 = Expiration date

 = Temperature limitation

 = Instruction for use

### Bibliografia

Rohle G, Siekmann L. Quality assurance of quantitative determination. In: Thomas L., editor Clinical laboratory diagnostics, 1st ed. Frankfurt: TH-Books Verlagsgesellschaft, 1998.p. 1393-1401.

REVISIONE	DATA	MODIFICA
Rev.A	10/2022	New issue for IVDR adjustment Regulation (EU) 2017/746

Scavo Diagnostics International  
Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)  
Phone +39 0577 39041 - Fax +39 0577 390 444

# Normal control serum - Instruction for use

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- Reattivi per sistemi analitici I.S.E. S.r.l.