 CUSTOMISED SOLUTIONS FOR YOUR LABORATORY	EN / IT	<b>HbA1c Control High</b>		<b>Instruction for Use</b> <b>Istruzioni per l'Uso</b>	 IVD	RevC - 15-11-22
		REF	R0400000004			

## Hemoglobin A1c Control High / HbA1c Controllo Alto

<b>LOT</b>	Expiry date Data di scadenza 
------------	---

### Intended Use

Accuracy control in the elevated range for quantitative determination of HbA1c in human blood by turbidimetric immunoassay.

**For professional in vitro diagnostic use only.**

### Composition

The HbA1c Control High is a hemolysate prepared from packed human erythrocytes, lyophilized and stabilised.

### Package

REF	R0400000004	1 x 0.5 mL
-----	-------------	------------

### Reagent Preparation

Accurately add 0.5 mL of distilled water to dissolve the contents. Gently mix for 30 minutes, or until all material has dissolved.

### Storage and Stability

The expiry date of the product at + 2 to + 8°C is listed on the label.

The reconstituted control is stable for 30 days if stored tightly closed at + 2 to + 8°C.

### Procedure

The lyophilized HbA1c control should be treated in the same manner as patient samples regarding the hemolysate procedure. Follow the directions that accompany the instrument and the reagent kit used in the assay for specific instrument calibration procedures.

### Precautions and Warnings

- In vitro diagnostic use only.
- Refer to the safety data sheets (SDS) and take the necessary precautions for the use with laboratory reagents.
- All blood donor units comprising the source material used in the manufacture of this product have been tested by approved FDA methods and found to be non-reactive for HBsAg, HCV and anti-HIV 1 & 2. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for bio hazardous materials since absence of infectious agents can never be proven.
- Do not use after expiry date.

### Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

### Reporting of serious incidents

The user must report (through the distributor) any serious accident occurring in relation to the device to both the manufacturer and the competent authority of the European Union Member State in which the user and / or patient is established. For other jurisdictions, reports of serious incidents must be produced in accordance with regulatory requirements.

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