 CUSTOMISED SOLUTIONS FOR YOUR LABORATORY	EN / IT	Anti-streptolysin (O) Standard High		Instruction for Use Istruzioni per l'Uso		 	Rev.C - 15-11-22
		REF	R1300001101				

Anti-streptolysin (O) Standard High Standard Alto Anti-streptolisina (O)

LOT	Expiry date Data di scadenza 
-----	-------------------------------------------------------------------------------------------------------------------

Intended Use

Preparation of reference curves for quantitative immunochemical determination of Anti-streptolysin O (ASL) in human serum.

For professional in vitro diagnostic use only.

Composition

The reagent is prepared by diluting human serum which contains high levels of ASLO with phosphate buffered saline containing 1 g % bovine serum albumin.

Preservative: 0.095 g % sodium azide.

Package

REF	R1300001101	Ready for Use	1 x 1 mL
-----	-------------	---------------	----------

Storage and Stability

The expiry date of the product at + 2 to + 8°C is listed on the label.

After first opening the container, the serum can be used for at least 6 weeks if stored tightly closed at + 2 to + 8°C after use.

Do not freeze.

Procedure

Refer to the IFU of the ASL(O) reagent used.

Precautions and warnings

- In vitro diagnostic use only.
- Refer to the safety data sheets (SDS) and take the necessary precautions for the use with laboratory reagents.
- All blood donor units comprising the source material used in the manufacture of this product have been tested by approved FDA methods and found to be non-reactive for HBsAg, HCV and anti-HIV 1 & 2. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for bio hazardous materials since absence of infectious agents can never be proven.
- Do not use after expiry date.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Reporting of serious incidents

The user must report (through the distributor) any serious accident occurring in relation to the device to both the manufacturer and the competent authority of the European Union Member State in which the user and / or patient is established. For other jurisdictions, reports of serious incidents must be produced in accordance with regulatory requirements.

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