 CUSTOMISED SOLUTIONS FOR YOUR LABORATORY	EN / IT	β2-Microglobulin Standard Set		Instruction for Use Istruzioni per l'Uso		 	RevC - 15-11-22
		REF	R1300002601				

β2-Microglobulin Standard Set / Calibratori β2-Microglobulina

LOT	Expiry date Data di scadenza 
-----	---

Intended Use

Preparation of reference curves for quantitative immunochemical determination of β2-Microglobulin (BMG).

For professional in vitro diagnostic use only.

Composition

BMG standard is a dilution of human serum spiked with purified β2-Microglobulin from human urine.

Preservative: 0.095 % sodium azide.

Package

REF	R1300002601		
Level 1, Level 2, Level 3, Level 4	Ready to Use	4x1 mL	

Storage and Stability

The expiry date of the product at + 2°C to + 8°C is listed on the label. After first opening the container, the serum can be used for at least 6 weeks if stored tightly closed at + 2°C to + 8°C after use.

Do not freeze.

Procedure

Refer to the IFU of the β2-Microglobulin reagent used.

Precautions and Warnings

- In vitro diagnostic use only.
- Refer to the safety data sheets (SDS) and take the necessary precautions for the use with laboratory reagents.
- Each individual donation intended for use in manufacture of protein control serum was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV1 and HIV 2 by FDA required test. Only donations with negative findings were used for its manufacture. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
- Do not use after expiry date.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Reporting of serious incidents

The user must report (through the distributor) any serious accident occurring in relation to the device to both the manufacturer and the competent authority of the European Union Member State in which the user and / or patient is established. For other jurisdictions, reports of serious incidents must be produced in accordance with regulatory requirements.

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