



<b>ISE</b>	<b>Transferrin</b>		<b>Instructions For Use</b>	 
	<b>REF R3330000060</b>		<b>100 Tests</b>	
	<b>R1</b> 1x 14.3mL	<b>R2</b> 1x 4.3mL	Doc Control: R3330000060_1.1A.E	

### Intended Use

Transferrin is a diagnostic assay for the quantitative detection of transferrin in human stool samples and is intended for use on the I.S.E. S.r.l. **Hemo One** analyser.

*For professional in vitro diagnostic use only.*

### Diagnostics Application

Transferrin is a blood-derived component that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases. Transferrin is stable in faeces and a good marker to detect loss of blood from the upper and lower intestine (gastrointestinal bleeding).

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme and heme-derived porphyrins in amounts that depend on the site, the amount of bleeding and the transit time through the gut.

### Method

Transferrin latex turbidimetric assay is based on agglutination reactions. These involve in vitro aggregation of microscopic latex particles. This aggregation consists in the specific reaction between antigen and antibodies, antigen contained in the sample and the antibodies anti-antigen coated on polystyrene latex particles. The sample is mixed with a suspension containing antibodies against the antigen bound to latex particles. If antigen is present in the sample, it will react with the antibodies and form an aggregate. If no antigen is present in the sample the mixture will keep its appearance as a smooth suspension. Such turbidity is measured as an increase in absorbance at the determinate wave and is proportional to the quantity of antigen contained in the sample.

### Reagents Provided

<b>R1</b>	<b>Reagent 1:</b>	1 Vials x 14.3 mL	Ready to Use
<b>R2</b>	<b>Reagent 2:</b>	1 vials x 4.3 mL	Ready to Use

### Stability and Storage

The reagents are stable until expiry date on the label when stored at 2-8°C. **DO NOT FREEZE**

After opening the reagents are stable for 1 month at 2-8°C.


### Reagents required but not supplied.

Reference	Description
R3330000061	Universal Stool Extraction Vials
R3330000067	Standard Set Transferrin
R3330000071	Control Set Transferrin

### Sample collection and preparation

Samples should be collected utilising the Universal Stool Extraction Vials. Homogenize stool samples as thoroughly as possible prior to preparation.

- R3330000061 Universal Stool Extraction Vials

 Consult the Universal Stool Extraction Vial instructions for use for correct stool sample extraction.

The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing.

### ASSAY PROCEDURE

#### Application parameters

The Application parameters are included in the Hemo One software. In the event of a missing the Transferrin protocol please contact your authorised ISE representative.

Any application not explicitly approved by I.S.E. S.r.l. cannot be guaranteed in terms of performance.

#### Hemo One Use

 Consult the I.S.E. S.r.l. Hemo One user manual for instructions on the proper use of the analyser.

#### Material Preparation

- Ensure samples are collected correctly in accordance with the sample preparation procedure.
- Allow reagents and stool samples to reach room temperature (15-30°C) prior to testing.
- Allow any calibrator or control material to reach room temperature (15-30°C) prior to testing.

#### Calibration curve establishment

A calibration curve must be performed prior to running samples. The curve must be validated with controls.

#### Calibration Stability

The calibration is recommended to be run every 7 days on the Hemo One analyser.

#### Quality control

For quality control purposes only use the following materials:

- R3330000071 Control Set Transferrin



Concentration is indicated on the label of the vial. The use of control materials at two different concentrations is recommended in order to verify test precision across the measuring range. Control frequency should be run in accordance with the laboratory's quality management system.

If the obtained results are out of the tolerance range;

- Ensure all materials are not expired
- Ensure all materials have been stored and prepared correctly
- Perform a calibration

In the event of controls results not meeting the defined tolerances please contact you local ISE authorised representative for support.



<b>ISE</b>	<b>Transferrin</b>		<b>Instructions For Use</b>	 
	<b>REF</b> R33300000060		<b>100 Tests</b>	
	<b>R1</b> 1x 14.3mL	<b>R2</b> 1x 4.3mL	Doc Control: R33300000060_1.1A.E	

## ANALYTICAL CHARACTERISTICS / PERFORMANCE

### Results

The results are evaluated automatically by the Hemo One and presented in ng hTf/mL.

### NORMAL VALUES

Positive results: A positive cut off should be determined by the clinical laboratory

Recommended cut off: 10 ng hTf/mL

### Method Comparison

Results obtained with the ISE Transferrin were compared with an immunochromatographic test (CerTest Transferrin, CerTest). The results were as follows:

	Sensitivity	Specificity
ISE Transferrin vs certest	<b>95%</b>	<b>&gt; 99%</b>

### Linearity

Using the calibrator kit the Transferrin is linear in the calibration range of **0-250 hTf/mL**.

### Limit of detection (LOD):

Limit of detection (LOD): **5 ng hTf/mL**

The lower limit of detection of Transferrin was determined on 20 samples and 2 sample replicates as the mean value + 2 SD

### Limit of quantification (LOQ):

Limit of quantification (LOQ): **7 ng hTf/mL**

The lower limit of quantification is defined as the lowest actual amount of analysis that can be reliably detected; imprecision is < 20% as CV%

### Precision

Precision was assessed with 3 different controls

Precision	Low	Medium	High
Number	20	20	20
Mean (ng/mL)	16.1	80.6	197.2
SD	1.2	3.6	10.3
CV (%)	7.4	4.7	5.2

### Prozone effect








Studies have been made up to a concentration of 2 µg of transferrin/mL and no false negative results have been observed. Studies using higher concentrations have not been carried out.

### Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

The product is in conformity with D.L: 8 September 2000, no. 332 "Actuation of the directive 98/79/EC regarding in vitro medical diagnostic devices".

### Symbols on labels and packaging

	In vitro diagnostic medical device
	Catalog Number
	Lot or batch number
	Manufacturer
	Expiry date
	Temperature limitation
	Consult Instructions for use

