

Diagnostic reagent for quantitative in vitro determination of gamma-glutamyltransferase (gamma-GT) in human serum or plasma on photometric systems

TEST PARAMETERS


Method: Colorimetric, kinetic, increasing reaction Szasz, standardized to IFCC
 Wavelength: 405 nm (400 – 420 nm)
 Temperature: 37 °C
 Sample: Serum, heparin plasma
 Linearity: up to 1200 U/L on automated systems
 Sensitivity: The lower limit of detection is 2 U/L

SUMMARY

Gamma-glutamyltransferase (gamma-GT / GGT), also called gamma-glutamyltranspeptidase, is an enzyme present in liver and bile duct which is the most sensitive indicator of hepatobiliary diseases. Because of a high negative predictive value for these diseases, the measurement of gamma-GT is widely used to rule out a hepatic or biliary origin. Together with other enzymes such as alanine aminotransferase (ALT), aspartate aminotransferase (AST) and cholinesterase gamma-GT is a valuable tool for the differential diagnosis in liver diseases. [1]

TEST PRINCIPLE

Kinetic photometric test according to Szasz/Persjin [2]. The test has also been standardized to the method according to IFCC [4]. Gamma-GT catalyzes the transfer of glutamic acid to acceptors like glycylglycine. This process releases 5-amino-2-nitrobenzoate, which can be measured at 405 nm. The increase in absorbance at this wavelength is directly related to the activity of gamma-GT.


REAGENT COMPOSITION

COMPONENTS	CONCENTRATION
Reagent 1:	
TRIS, pH 8,28	135 mmol/L
Glycylglycine	135 mmol/L
Reagent 2:	
L-Gamma-glutamyl-3-carboxy-4-nitroanilide, pH 6,00	22 mmol/L

REAGENT PREPARATION
Substrate Start:

The reagents are ready for use.

Sample Start:

Mix 4 parts of Reagent 1 + 1 part of Reagent 2

(= working reagent).

REAGENT STABILITY AND STORAGE

Conditions: Protect from light
 Close immediately after use. Avoid contamination.
 Do not freeze!

Substrate Start:

Storage: at 2 – 8 °C

Stability: up to the expiration date indicated on labels

Sample Start (working reagent):

Storage: at 2 – 8 °C at 15 – 25 °C

Stability: 4 weeks 5 days

The working reagent must be protected from light.

SAMPLE STABILITY AND STORAGE

Serum, heparin plasma:

Stability [6]: at -20 to +25°C at least 1 week

Freeze only once!

Discard contaminated specimens.

MATERIALS REQUIRED BUT NOT PROVIDED

NaCl solution (9 g/L)

General laboratory equipment

MANUAL TEST PROCEDURE

Bring reagents and samples to room temperature.

Substrate start 37 °C

Pipette into test tubes:	Blank	Sample
Reagent 1	1000 µL	1000 µL
Distilled water	100 µL	
Sample		100 µL
Mix. Incubate for approximately 1 minute. Then add:		
Reagent 2	250 µL	250 µL
Mix. Read initial absorbance after 1 minute and start a timer. Read absorbance again after exactly 1, 2 and 3 minutes. Determine ΔA/min. during the linear part of the assay. ΔA/min = [ΔA/min sample] – [ΔA/min blank]		

Sample start 37 °C

Pipette into test tubes	Blank	Sample
Working reagent	1000 µL	1000 µL
Distilled water	100 µL	
Sample		100 µL
Mix. Read initial absorbance after 1 minute and start a timer. Read absorbance again after exactly 1, 2 and 3 minutes. Determine ΔA/min. during the linear part of the assay. ΔA/min = [ΔA/min sample] – [ΔA/min blank]		

CALCULATION

With factor: (light path 1 cm)

gamma-GT (U/L) = ΔA/min x factor

Factors (37°C):

	Szasz	IFCC
For Substrate start:	1421	1606
For Sample start:	1158	1309

With calibrator:

$$\text{GGT [U/L]} = \frac{\Delta\text{A/min Sample}}{\Delta\text{A/min Calibrator}} \times \text{activity calibrator [U/L]}$$

Results according to Szasz or IFCC resp. are obtained by use of the calibrator value given for the Szasz method, or for the IFCC method respectively.

UNIT CONVERSION

U/L x 0.01667 = µkatal/L



REFERENCE RANGE *

According to Szasz [5]:

	U/L	μkat/L
Women	< 32	0.53
Men	< 49	0.82

According to IFCC:

	Female [U/L]	Male [U/L]
Adults [4]	< 38	< 55
Children/adolescents: [1]		
1 day – 6 months	15 – 132	12 – 122
6 months – 1 year	1 – 39	1 – 39
1 – 12 years	4 – 22	3 – 22
13 – 18 years	4 – 24	2 – 42

	Female [μkat/L]	Male [μkat/L]
Adults [4]	< 0.63	< 0.92
Children/adolescents: [1]		
1 day – 6 months	0.250 – 2.20	0.200 – 2.03
6 months – 1 year	0.017 – 0.651	0.017 – 0.651
1 – 12 years	0.067 – 0.367	0.050 – 0.367
13 – 18 years	0.067 – 0.401	0.033 – 0.701

* Each laboratory should check if the reference ranges are transferable to its own patient population and determine own reference ranges if necessary.

PERFORMANCE CHARACTERISTICS

LINEARITY, MEASURING RANGE

On automated systems, the test is suitable for the determination of gamma-GT activities up to 1200 U/L.

In case of a manual procedure, the assay is suitable for gamma-GT activities which correspond to a maximum of $\Delta A/\text{min} = 0.20$.

If such values are exceeded, the samples should be diluted 1+5 with NaCl solution (9 g/L) and the results multiplied by 6.

SENSITIVITY/LIMIT OF DETECTION

The lower limit of detection is 2 U/L

PRECISION

Intra-assay n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	39.9	0.99	2.48
Sample 2	73.6	0.85	1.16
Sample 3	206	1.32	0.64

Inter-assay n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	41.5	0.62	1.49
Sample 2	72.3	0.61	0.85
Sample 3	204	0.74	0.36

SPECIFICITY/INTERFERENCES

No interference up to:

Ascorbic acid	30 mg/dL
Bilirubin	40 mg/dL
Hemoglobin	400 mg/dL
Triglycerides	2000 mg/dL

For further information on interfering substances refer to Young DS [7].

METHOD COMPARISON

A comparison of our gamma-GT, standardized to IFCC (y) with the IFCC reference reagent (x) using 51 samples gave following results:

$$y = 1.005x - 0.741 \text{ U/L}; r = 0.999$$

A comparison of our Gamma-GT, Szasz (y) with a commercially available test according to Szasz (x) using 51 samples gave following results: $y = 0.996x + 1.354 \text{ U/L}; r = 1.000$.

CALIBRATION

The use of a gamma-GT calibrator is optional.

We recommend the our multi calibration serum. In case this calibrator is used, use the according calibrator value for the Szasz method respectively for the IFCC method.

For calculation according to IFCC, standardization was performed against the original IFCC formulation.

QUALITY CONTROL

Control sera with gamma-GT values determined by this method can be used.

We recommend the our serum controls (control serum with values in the normal and abnormal range)

Each laboratory should establish corrective action in case of deviations in control recovery.

WARNINGS AND PRECAUTIONS

- The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
- In very rare cases, samples of patients with gammopathy might give falsified results [8].
- Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents.
- For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examinations and other findings.
- For professional use only!

WASTE MANAGEMENT

Please refer to local legal requirements.

REFERENCES

- Thomas L. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft;1998.p.80-6.
- Persijn JP, van der Silk W. A new method for the determination of gamma-glutamyltransferase in serum. J Clin Chem Clin Biochem 1976; 14:421-7.
- Szasz G. Gamma-Glutamyltranspeptidase. In: Bergmeyer HU. Methoden der enzymatischen Analyse. Weinheim: Verlag Chemie, 1974. p. 757.
- Schumann G, Bonora R, Ceriotti F, Férard G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37°C. Part 5: Reference procedure for the measurement of catalytic concentration of gamma-glutamyltransferase. Clin Chem Lab Med 2002; 40:734-8
- Fischbach F, Zawta B. Age-dependent reference limits of several enzymes in plasma at different measuring temperatures. Klin Lab 1992;38:555-61.
- Guder WG, Zawta B et al. The quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 30-1.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
- Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007; 45(9): 1240-1243.

Symbols on labels and packaging

= In vitro diagnostic medical device

= Catalog Number

= Lot Number

= Manufacturer

= Expiration date

= Temperature limitation

= Instruction for use

