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ISE S.r.l.

CUSTOMISED SOLUTIONS
FOR YOUR LABORATORY

4 x 3 mL •  R0300000006

Intended use

For calibration of the I.S.E. test methods on all the I.S.E. Clinical Chemistry Analyzers. The reagents and methods described are manufactured and controlled by Scavo Diagnostics International. For professional use only.

Composition

The Calibration Serum is obtained from defibrinated, freeze-dried human plasma. The declared nature of the product denotes the potential risk involved and it must therefore be used, adopting the same precautions used in handling the test samples.

Reconstitution of the lyophilized serum

Remove the cap with care, to avoid losing traces of the freeze-dried product. Add exactly 3 mL of bidistilled water, recap, mix by inversion a few times and leave to rest for about 30 minutes. Before use, mix again by inversion, avoiding foam formation.

Important: when calibrating for Alkaline Phosphatase (ALP), wait a further 60 minutes until the enzyme is completely reactivated.

Storage and stability of the lyophilized serum

Store at 2 - 8°C. The stability of the serum before reconstitution is reported on the vial label. Stability of the reconstituted serum.

at 18 - 25°C 8 hours
at 2 - 8°C 2 days
at (-15) – (-25)°C 1 month (If thawed only once)

Origin of the biological components

Analyte	Origin
ALT (GPT)	Pig heart
AST (GOT)	Pig heart
Albumin	Human serum
Amylase	Human saliva/pig pancreas
Cholesterol	Bovine plasma
Cholinesterase	Human serum
CK	Pig plasma
Y-GT	Pig kidney
LDH	Pig heart
ALP (Alkaline phosphatase)	Calf intestine
ACP (acid phosphatase)	Human prostate/potato
Total Proteins	Human serum

Precautions and warnings

The sera have been tested and proved negative for the presence of HbsAg, and also for anti-HIV 1 and 2 and anti-HCV antibodies. However, the material must be considered potentially dangerous and handled in the same way as the test samples.

Reporting of serious incidents

Please inform the manufacturer (through your distributor) and the competent authority of the member state of the European Union in which the user and / or patient is established, of the cases of serious accident occurring in relation to the device. For other jurisdictions, reporting of serious incidents must be made in accordance with regulatory requirements. By reporting serious incidents, you help provide more information about the safety of the in vitro diagnostic medical device.

Assigned values and traceability

The theoretical values have been calculated on the basis of the results obtained on I.S.E. biochemical analyzers using I.S.E. reagents. The results of the concentrations are obtained with reference to International Reference Standards. The traceability of the values, with reference to each diagnostic method, is reported on the relative product information sheet. The use of procedures or reagents other than those indicated may lead to differences between the expected values and those obtained.

Materials included in the kit

- Freeze-dried serum as described
- Product information sheet reporting the theoretical values for the calibrator

Materials required but not included in the kit

- Distilled or deionized water
- Reagents for I.S.E. analytical systems

Disposal

Disposal must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Manufacturer:

Scavo Diagnostics International
Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)
Phone +39 0577 39041 - Fax +39 0577 390 444

Distributor:

I.S.E S.r.l.
Via Delle Driadi, 45 – 00133 Roma
Tel. +39 077 4579365; FAX +39 077 4579305
E-mail: info@logotech-ise.com
www.logotech-ise.com

Symbols used on labels and packaging

 = In vitro diagnostic medical device

 = Catalogue Number

 = Lot Number

 = Manufacturer

 = Expiration date

 = Storage temperature

 = Instruction for use

REVISIONE	DATA	MODIFICA
Rev.A	10/2022	New issue for IVDR adjustment Regulation (EU) 2017/746

Scavo Diagnostics International
Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)
Phone +39 0577 39041 - Fax +39 0577 390 444

Multicalibrator – Instruction for use

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