 CUSTOMISED SOLUTIONS FOR YOUR LABORATORY	EN / IT	HbA1c Calibrator (5 levels)		Instruction for Use Istruzioni per l'Uso		 IVD	Rev.C - 15-11-22
		REF	R0300000003				

Hemoglobin A1c Calibrator / Calibratore HbA1c (5 livelli)

LOT	Expiry date Data di scadenza 
-----	---

Intended Use

Preparation of reference curves for the quantitative determination of HbA1c in human blood by turbidimetric immunoassay.

For professional in vitro diagnostic use only.

Composition

The HbA1c Standard Set Level 1 is a liquid buffer with preservative: 0.095 % sodium azide. The HbA1c Standard Set Level 2, Level 3, Level 4 and Level 5 are hemolysates prepared from packed human erythrocytes, lyophilized and stabilised.

Package

REF	R0300000003	
Level 1	Ready to Use	2 mL
Level 2, Level 3, Level 4, Level 5	To be reconstituted	0.5 mL

Reagent Preparation

HbA1c Standard Set Level 1: ready for use.

HbA1c Standard Set Level 2, Level 3, Level 4 and Level 5: accurately add 0.5 mL of distilled water to dissolve the contents. Gently mix for 10 minutes, or until all material has dissolved.

Storage and Stability

The expiry date of the product at + 2 to + 8°C is listed on the label.

The reconstituted standard set is stable for 30 days if stored tightly closed at + 2 to + 8°C.

Procedure

The lyophilized HbA1c standard set will produce a calibration curve that will be stable for at least 2 weeks on most of analysers. The standards should be treated in the same manner as patient samples regarding the hemolysate procedure. Follow the directions that accompany the instrument and the reagent kit used in the assay for specific instrument calibration procedures.

Use fresh standards to prepare HbA1c hemolysate standards 2 to 5:

Lysis (for A-R0000000030)

Dispense 2 mL Hemolysis Solution for HbA1c into a test tube. Place 20 µL of well mixed calibrator. Mix. Allow to rest for 5 minutes or until complete lysis is evident.

Precautions and Warnings

- In vitro diagnostic use only.
- Refer to the safety data sheets (SDS) and take the necessary precautions for the use with laboratory reagents.
- All blood donor units comprising the source material used in the manufacture of this product have been tested by approved FDA methods and found to be non-reactive for HBsAg, HCV and anti-HIV 1 & 2. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for bio hazardous materials since absence of infectious agents can never be proven.
- Do not use after expiry date.

Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Reporting of serious incidents

The user must report (through the distributor) any serious accident occurring in relation to the device to both the manufacturer and the competent authority of the European Union Member State in which the user and / or patient is established. For other jurisdictions, reports of serious incidents must be produced in accordance with regulatory requirements.

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