ISE S.r.l. CUSTOMISED SOLUTIONS FOR YOUR LABORATORY	EN / IT	Ferritin Standard Set		Instruction for Use Istruzioni per l'Uso	 IVD	Rev.C - 15-11-22
		REF	R1300002901			

Ferritin Standard Set / Calibratori Ferritina

LOT	Expiry date Data di scadenza		
------------	---------------------------------	--	--

Intended Use

Preparation of reference curves for quantitative determination of Ferritin in human serum by turbidimetric immunoassay. The standard set has to be used with the Ferritin KIT 2nd Generation.

For professional in vitro diagnostic use only.

Composition

Ferritin Standard Set are dilutions of human plasma, liquid stabilised and filtered through 0.2 µ.

Preservative: Sodium azide (0.95 g/L).

Package

REF	R1300002901		
Level 1, Level 2, Level 3, Level 4, Level 5	Ready to Use	5x1 mL	

Storage and Stability

The expiry date of the product at + 2°C to + 8°C is listed on the label.

The reconstituted control is stable for 6 weeks if stored tightly closed at + 2°C to + 8°C.

Do not freeze.

Procedure

Refer to the IFU of the Ferritin reagent used.

Precautions and Warnings

- In vitro diagnostic use only.
- Refer to the safety data sheets (SDS) and take the necessary precautions for the use with laboratory reagents.
- All blood donor units comprising the source material used in the manufacture of this product have been tested by approved FDA methods and found to be non-reactive for HBsAg, HCV and anti-HIV 1 & 2. Nevertheless every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for bio hazardous materials since absence of infectious agents can never be proven.
- Do not use after expiry date.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

Disposal of reagent

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

Reporting of serious incidents

The user must report (through the distributor) any serious accident occurring in relation to the device to both the manufacturer and the competent authority of the European Union Member State in which the user and / or patient is established. For other jurisdictions, reports of serious incidents must be produced in accordance with regulatory requirements.

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