

R1: 4 x 46 mL – R2: 4 x 12 mL •  A-R0100000201R1: 3 x 12,6 mL – R2: 1 x 9,5 mL •  R3330000010**INTENDED USE**

Product for use in the quantitative determination in vitro of the concentration of the Direct Bilirubin in human serum or plasma. The results of the test must always be interpreted in conjunction with the clinical context. For professional use only.

**CLINICAL SIGNIFICANCE**

Approximately 80-85 % of the bilirubin produced is derived from the heme moiety of the haemoglobin released from aging erythrocytes in the reticuloendothelial cells. Bilirubin, bound to albumin, is transported into the liver where it is rapidly conjugated with glucuronide to increase its solubility. Then it is excreted into biliary canaliculi, and hydrolyzed in the gastrointestinal tract.

Unconjugated bilirubin serum concentration increases in case of overproduction of bilirubin (acute or chronic haemolytic anemias) and in case of disorders of bilirubin metabolism and transport defects (impaired uptake by liver cells: Gilbert's syndrome; defects in the conjugation reaction: Crigler-Najjar syndrome). Reduced excretion (hepatocellular damage: hepatitis, cirrhosis, Dubin-Johnson and Rotor syndrome) and obstruction to the flow of bile (most often produced by gallstones or by tumors) induce an important elevation of conjugated bilirubin and in a minor extent an increase of unconjugated bilirubin (conjugated hyperbilirubinemia).

**PRINCIPLE**

Mally – Evelyn modified. End point. Sulphanilic acid reacts with sodium nitrite to form diazotized sulphanilic acid. In the presence of accelerator (centrimide), conjugated and unconjugated bilirubin react with diazotized sulphanilic acid to form azobilirubin (Bilirubin Total 4+1). In the absence of accelerator, only conjugated bilirubin reacts (Bilirubin direct 4+1). The increase of absorbance at 550 nm is proportional to bilirubin concentration.

Sulphanilic acid + NaNO<sub>2</sub> → Diazotized sulphanilic acid

Bilirubine + Diazotized sulphanilic acid → Azobilirubin

**REAGENTS**

A-R0100000201 – R1: 4 x 46 mL – R2: 4 x 12 mL

**Reagent 1:** n° 4 vials x 46,0 mL ready to use

**Reagent 2:** n° 4 vials x 12,0 mL ready to use

R3330000010 – R1: 3 x 12,6 mL – R2: 1 x 9,5 mL

**Reagent 1:** n° 3 vials x 12,6 mL ready to use

**Reagent 2:** n° 1 vial x 9,5 mL ready to use

**Concentrations**

Reagent 1:		
	Conc.	U.M.
Sulphanilic Acid	29.0	mmol/L
Reagent 2:		
Sodium Nitrite	11.0	mmol/L

**Precautions**

**Reagent 1 contains sulphanilic acid.** May produce an allergic reaction. Consult the Safety Data Sheet. In addition to risk indications related to the active components, reagents may contain inactive components such as preservatives and detergents. The total concentration of these components is lower than the limits reported in the EC 1272/2008 regulation and subsequent amendments and additions. However, It is recommended to handle reagents according to the standards of good laboratory practice.

**Reports of serious incidents**

In the event of a serious incident in relation to the use of the device, please inform the manufacturer (via your distributor) and the competent authority of the European Union member state where the incident occurred. For other jurisdictions, reporting must be made in accordance with regulatory requirements. Reporting serious incidents helps provide more information about the safety of the diagnostic medical device.

**Storage and stability**

Store at 2-8°C and protect from direct light. When correctly stored, the reagents are stable up the expiry date reported on the label. A slight variation in the colour of the reagent may occur between batches, but this has no effect on the test results. After opening, the vial R1 and R2 are stable 28 days if recapped immediately and protected from contamination, evaporation, direct light and stored at correct temperature

**SAMPLE COLLECTION****Type of sample and storage**

Fresh, non-haemolyzed serum or plasma samples should be used. Store away from the light. The samples must be tested within 2 hours if stored at room temperature or within

12 hours if stored at 2 - 8°C. If frozen between - 15°C and - 20°C the samples are stable for 3 - 4 months (4).

**Precaution**

All human samples must be handled and disposed of as potentially infectious materials.

**PROCEDURE****Quality control**

Human control serum with known levels of Direct Bilirubin is commercially available for quality control purposes. Data sheets are included, listing the values and the confidence limits. Normal and abnormal control sera are available as "Normal control serum" code R0400000006 and "Pathological control serum" code R0400000106. The values obtained must be contained within the acceptability range. In case of incorrect results check the following points:

Cleanliness of glassware.

- Wavelength setting.
- Expiration date of reagents.

**Automation**

Although this device has been developed and manufactured to be used by manual method and instrumental systems, it can also be used in combination with other instrumental devices that are able to meet the requirements indicated in the paragraph "Reaction conditions / Technique".

All applications that are not explicitly approved cannot be guaranteed in terms of performance and will therefore have to be evaluated by the user

**Calibration**

For calibration use the "Multicalibrator code R0300000006. For the use on ISE srl instruments Calibration every 10 days is recommended.

**Traceability:**

The Direct Bilirubin value is reported in the package insert supplied with the "Multicalibrator".

**Reagents included in the kit**

The reagent is described above.

**Materials required but not supplied in the kit**

Calibrators and controls.

**Reaction conditions**

Wavelength (primary): 550 nm

Temperature: 37°C

**Technique**

Bring the reagents to the reaction temperatures

	U.M.	Calibration Serum	Sample	Blank
<b>Reagent A</b>	µL	800	800	800
<b>Calibration Serum</b>	µL	100	-	-
<b>Sample</b>	µL	-	100	-
<b>Blank</b>	µL	-	-	100
Mix gently and incubate at 37°C for 5 min. Read absorbance A1 at 550 and 700 nm, then add. DA1= A1 550nm - A1 700nm				
<b>Reagent B</b>	µL	200	200	200
Mix gently and incubate at 37°C for 50 sec. Read absorbance A2 at 550 and 700 nm, then add. DA2= A2 550nm - A2 700nm				

**Reading**

Read the absorbances at 550 nm and 700 nm of the sample and the calibrator, subtracting the absorbance of the reagent white. The final color is stable 1 hour away from direct light.

**Results**

The concentration of Direct Bilirubin is obtained through the following formula

$[A\Delta 2 - A\Delta 1 \times Fdil]_{sample} \times n$

$[A\Delta 2 - A\Delta 1 \times Fdil]_{calibrator}$

Fdil = Dilution factor (0.81)

N = Calibrator concentration

Conversion factor mg/ dl x 17.1 = µmol/L



R1: 4 x 46 mL – R2: 4 x 12 mL •  A-R0100000201R1: 3 x 12,6 mL – R2: 1 x 9,5 mL •  R3330000010**NORMAL VALUES**

Serum or plasma:

Up to 0.2 mg/dL (3.4 µmol/L).

Each laboratory must establish its own normal values on the basis of its local population.

**ANALYTICAL CHARACTERISTICS / PERFORMANCE****Linearity**

The measurement range is 0.08 to 10.55 mg/dL (1,4 – 180,4 µmol/L) of Direct Bilirubin in the sample. For higher concentrations, repeat the test on the sample diluted 1:5 in physiological solution then multiply the result by the dilution factor. This procedure extends the measuring range up to 50 mg/dL (855,2µmol/L)

**Interference**

No significant interference was detected in the presence of triglycerides (up to 2000 mg/dl), Hemoglobin (up to 125 mg/dl); Acetaminophene (up to 30 mg/dl); Ascorbic acid (up to 0.5 mg/dl); Acetyl salicylic acid (up to 200 mg/dl). Do not use hemolyzed samples.

**Accuracy-Recovery**

Accuracy studies were conducted on normal samples to which amounts of pure bilirubin notes in the range 1 - 10 mg/dL were added. The data indicate a recovery of 99.95%.

**Precision of the method**

Within-run precision				
Range	U.M.	Mean	C.V. (%)	No.
Low	mg/dL	0.36	3.8	80
Medium	mg/dL	1.51	1.9	80
High	mg/dL	3.99	0.9	80
Between-run precision				
Range	U.M.	Mean	C.V. (%)	No.
Low	mg/dL	0.36	5.2	80
Medium	mg/dL	1.51	5.3	80
High	mg/dL	3.99	4.7	80

**Sensitivity**At  $\lambda$  546 nm a concentration of about 0.08 mg/dL of Direct Bilirubin in the conditions established for this test.**Comparative method**The method was compared with a similar method on the market. Samples analysed n° = 100 (with values between 0.09 and 10.52 mg/dL); Regression line  $y = 0.926x - 0.03$ ;  $r = 0.998$ .**Disposal of reagent**

Disposal of reagents must be performed in accordance with the EC regulations regarding waste, or the local national or regional legislation.

**Manufacturer:****Sclavo Diagnostics International**

Via Po 26-28 – Loc. Pian dei Mori – 53018 Sovicille (SI) (Italy)

Phone +39 0577 39041 - Fax +39 0577 390 444

**Distributor:****I.S.E S.r.l.**








Via Delle Driadi, 45 – 00133 Roma

Tel.+39 077 4579365; FAX +39 077 4579305

E-mail: info@logotech-ise.com

www.logotech-ise.com

**Symbols used in IFU and Packaging**

 In vitro diagnostic medical device vitro	 Manufacturer
 Catalogue Number	 Instruction for use
 Lot Number	 Temperature limitation
 Expiration date	

**References**

- Malloy HT, Evelyn KA. J Biol Chem 1937; 119: 481-490.
- Colombo JP, Peheim E, Kyburz S, Hoffmann JP. Clin Chim Acta 1974; 51: 217-219.
- Winsten S, Cehelyk B. Clin Chim Acta 1969; 25: 441-446.
- Martinek RG. Clin Chim Acta 1966; 13: 161-170.
- Annino JS. Clinical Chemistry: Principles and Procedures. Little, Brown and Company, Boston-Toronto 1960; 198-203.
- Doumas BT, PerryBW, Sasse EA, Straumfjord JV. Clin Chem 1973; 19:984-993.

REVISION	DATE	CHANGE
Rev.C	01/2024	Edit header and "Reagents" paragraph.

