



Hema One^{Plus}


Compact & Versatile
Closed Fully Automated Solution
for testing Metabolic &
Gastro-Intestinal (Fecal) Disorders

- Reliable • User-friendly • Fast

The specific clinical panels are designed to give you fast, direct and reliable information about the status of metabolic functionality

EMPOWER METABOLIC TESTING

- GASTRO-INTESTINAL (FECAL)
- DIABETES
- CARDIAC PANEL
- VITAMIN DEFICIENCIES
- KIDNEY FUNCTION
- METABOLIC DISORDERS




Features:

- Fecal extraction system
- 9 samples controls & calibrators positions
- Cartridge based solution - 16 cartridge positions
- HbA1c lysing on board
- Up to 120 tph for single-reagent test
- Up to 60 tph for bi-test
- Sample tubes: 12-13 mm Ø and 75-85 mm height, 3 mL sample cups
- 72 disposable cuvettes on board
- Integrated operator console

Hemo One^{Plus}

Compact & Versatile
Closed Fully Automated Solution
for testing Metabolic &
Gastro-Intestinal (Fecal) Disorders



11154-02-A



CUSTOMISED SOLUTIONS
FOR YOUR LABORATORY



To find more about our Hemo One Plus, Visit Us at:

www.logotech-ise.com

I.S.E. S.r.l. - Via Luigi Einaudi, 7, Guidonia Montecelio, 00012 Rome, Italy



REAGENT MENU

GASTRO-INTESTINAL (FECAL) PANEL

Our GI (fecal) panel allows for rapid and accurate automated testing in stool samples for the assessment of intestinal inflammation and infections.

HAEMOGLOBIN (FOB)
CALPROTECTIN

H. PYLORI
TRANSFERRIN

DIABETIC PANEL

Our diabetic panel provide rapid, accurate and comprehensive testing to diagnose the major diabetic parameters.

HbA1c (Whole blood)

MICRO-ALBUMIN

CARDIAC PANEL

Our panel checks immediately several cardiac biomarkers. It helps to determinate a person's risk of cardiac ischemia, insufficient blood flow to the heart, certain thrombosis events and more...

D-DIMER
CK-MB
CK-NAK

LACTATE DEHYDROGENASE
hs- CRP

KIDNEY FUNCTION PANEL

Our panel is used as part of screening or assessment of risk to develop a kidney disease, or to monitor a patient with a known kidney symptoms (e.g. high blood pressure, diabetes, proteinuria). Panel can be used for monitoring the response to a treatment as well.

CYSTATIN-C

SODIUM
POTASSIUM

ADDITIONALITY AVAILABLE TESTS

VITAMIN D
VITAMIN B12 *
VITAMIN 9 *

* Coming soon

Customer Care Center:
info@logotech-ise.com

Tel: +39 0774 579 326
Fax: +39 0774 579 305

Technical Service:
service@logotech-ise.com

Tel: +39 0774 579 364
Fax: +39 0774 579 085

Sales and Marketing office:
sales@logotech-ise.com

Tel: +39 0774 579 334
Fax: +39 0774 579 305

Order handling and Logistics:
order@logotech-ise.com

Tel: +39 0774 579 365
Fax: +39 0774 579 085